

## **Product Specifications**

<b>Product number:</b>	SMA-3010
<b>Description:</b>	Smith Antigen Coated Microtiter Plate
<b>Amount:</b>	1X 96 wells
<b>Manufacturing method:</b>	Smith Antigen is absorbed onto each microtiter well at 1unit/well. The wells are post-coated with a stabilizing solution.
<b>Storage conditions:</b>	Store at 4° C in its re-sealable zip lock pouch with desiccant
<b>Expiration:</b>	See Label
<b>Instructions for use:</b>	For use with standard solid phase immunoassay procedures.
<b>Warnings:</b>	For research or further manufacturing use only.

QC0083  
SMA3010.DOC  
Ver 1.3  
Effective Date 9-5-2003