

Product Specifications

<i>Product number:</i>	<i>PAG-3010</i>
<i>Description:</i>	<i>Ribosomal P Antigen Coated Microtiter Plate</i>
<i>Amount:</i>	<i>1X 96 wells</i>
<i>Manufacturing method:</i>	<i>Ribosomal P Antigen is absorbed onto each microtiter well at 1unit/well. The wells are post-coated with a stabilizing solution.</i>
<i>Storage conditions:</i>	<i>Store at 4° C in its re-sealable zip lock pouch with desiccant</i>
<i>Expiration:</i>	<i>See Label</i>
<i>Instructions for use:</i>	<i>For use with standard solid phase immunoassay procedures.</i>
<i>Warnings:</i>	<i>For research or further manufacturing use only.</i>